

ActivForce™ Features

Benefits

• Universal design (right or left leg and either compartment)

• Eight Sizes (Circumference of the Knee):

12" - 13" XS	17" - 19" XL
13" - 14" S	19" - 21" XXL
14" - 15" M	21" - 22" 3XL
15" - 17" L	22" - 23" 4XL

• Light, low profile design

• Patented Townsend Motion TM6 Hinge

• Lightweight, breathable, non-neoprene sleeve with lycra posterior panel and wide 2-way thigh cuff opening

• Sleeve incorporates anti-migration silicone beads

• Three condylar pad thicknesses to choose from (3/16", 1/4", & 3/8")

• Allows the hinge to be inserted into the medial or lateral sleeve pocket for either leg

• Fits patients with varying leg sizes and shapes

• Can be worn under clothing for improved compliance

• Tracks anatomically with the knee for more natural movement than traditional polycentric hinge designs

• Meets L-1843 code requirements for a single hinge brace with adjustable extension/flexion control

• Easy for patient to put on brace and achieve compression to reduce swelling

• Stabilizes knee

• Keeps brace comfortably in place during periods of activity

• Allows patient to adjust and optimize corrective force